




ATM Qness GmbH
Emil-Reinert-Str. 2
57636 Mammelzen
Germany

Phone +49 2681 9539-0
E-Mail info@qatm.com

ATM Qness GmbH
Reitbauernweg 26
5440 Golling
Austria

Phone +43 6244 34393
E-Mail info.at@qatm.com

www.qatm.com

ELECTROLYTE POLISHING AND ETCHING OF A CAST ALUMINUM ALLOY

Electrolytic polishing and etching is an appropriate preparation method with less deformation for non-ferrous materials and materials with a single-phase microstructure that must later be analyzed under an EBSD detector.

Electrolytic polishing directly on the finished part can save time and costs by avoiding the normal grinding/polishing process.

A well-known application of electrolytic etching is the use of Barker solution to enable color etching in Aluminum alloys. The etching affects the individual grains of the microstructure depending on their orientation. If the etched samples are viewed under linearly polarized light, the differently oriented grains appear in various colors. With this method, the direction of growth in the investigated plane can be easily recognized and the size of the grains can be measured.

Objective:

The Aluminum sample should be prepared to show the dendritic structure and phase distribution. The preparation should include the following steps: electrolytic polishing and etching with Qetch 1000. Polishing was performed with the K1 electrolyte and etching with the Barker etchant.

Aluminum sample



ELECTROLYTIC POLISHING AND ETCHING




Electrolytic polisher and etcher
Qetch 1000

Device		Qetch 1000				
Step		ELECTROLYTE	VOLTAGE [V]	PUMPING POWER [%]	MASK [CM ²]	DURATION [MIN]
ELECTROLYTIC POLISHING						
	Polishing (electrolytic)	Electrolyte K1	20	55	2.5	0:30
ELECTROLYTIC ETCHING						
	Etching (electrolytic)	Barker	20	65	2.5	1:30
Notes						


The sample is placed directly into the QETCH 1000.



The structure after electrolytic etching with
the Barker etchant – 25 x



The dendritic microstructure of
an aluminum cast alloy – 25 x



Find out more at
www.qatm.com