MICRO DUROMETRI VICKERS / KNOOP / BRINELL / ROCKWELL

QNESS 60 A+ EVO

qatm.com/it/prodotti/test-di-durezza/micro-hardness-tester/qness-60-a-evo/

QNESS 60 A+ EVO

MICRO DUROMETRI

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The innovative „EVO“ model of this Vickers / Knoop / Brinell / Rockwell hardness tester promises perfect automation and comes with a high-precision positioning control for all three axis systems with glass scales.

Flexible XYZ progressions for thousands of analysis points allow for exceptionally high throughput without operator intervention.

QATM micro hardness testers seamlessly integrate with the Verder Scientific IoT platform, featuring remote real-time monitoring, live notifications, effortless backups and automatic software updates.

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
QNESS 60 A+ EVO

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The QATM-developed, in-house manufactured lens system sets new standards. As well as providing crystal clear image quality for hardness testing, Koehler illumination uses white LED light and motor-operated aperture shuttering to produce ideal contrast, even for high magnification images.

Experienced metallurgists agree the image quality provided by the Qness 10 / 60 A+ is comparable in all aspects with that of established sophisticated microscopes.

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QATM hardness testers accurately analyze according to all standard test methods and cover a wide spectrum. The electronically controlled, fully automated test cycles ensure fast, precise hardness testing, as well as fast method changes and automatic detection of the focal plane. With the Ready to Test package all Vickers test methods are possible with the standard scope of delivery.

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


NEW ROCKWELL TEST METHOD

In addition to the Vickers, Knoop and Brinell optical methods, Rockwell measurements can now also be carried out with the new micro hardness testers. A special Rockwell test module has been developed for this purpose.

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



In a few seconds the image of the sample is shot with the additional camera (field of view 52 x 39 mm). The image provides excellent navigational support within the software, particularly in combination with DOUBLE-VIEW TECHNOLOGY.

HIGH-RESOLUTION SAMPLE IMAGE (HRI)

If high-quality images of larger areas are required (e.g. for weld seam measurements), the area can be scanned using the HRI function. The Qpix Control 2 software automatically combines the individual images into one large overall image.

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All 3 axes are equipped with incremental measuring system. The axes and turret can be positioned to an accuracy of 4 μm , so the perfect solution for classic micro hardness testing. Series measurements can be repeatedly and accurately approached.

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The optional large slide doubles the test table surface area to 300 x 120 mm. If required, the test height can also be extended from 150 mm to 260 mm.

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Optimized testing parameters and shorter intervals for serial autofocus, brightness regulation and image analysis, facilitate unbeatable cycle times during everyday operation involving hardness testing devices of the new micro hardness tester product line; and it's even faster than the previous model.

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

NON-EMBEDDED SAMPLES

Components of almost all geometrical shapes can be fitted into the universal sample holder. Four clamping bolts can be set variably in various T-slots.

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
PRISMS

QATM prisms also enable round components to be tested with our devices. Benefit: Integration of the 3D model in the software automatically determines the center of the component and the highest point of each piece.


VICES

The clearly-organized, large and robust test room facilitates universality. Additionally, clamping samples straight into conventional vices reduces the effort required to prepare the sample, and expands the range of possible future test applications.

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


1 | LOAD SAMPLES

The machine moves automatically to sample holder height. Image of sample is taken automatically.

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
SRS - Speedy Row Set-up: Drag the row of test points to the desired position. The serial Auto-Snap function corrects the starting point of the test row automatically.

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
INNOVATIVE SOFTWARE FUNCTIONS


FULLY AUTOMATIC HARDNESS TESTING

Several progressions and samples are created and completed "unmanned" (i.e. 60 progressions on 8 different samples in one test run).


AUTOMATIC MINIMUM DISTANCE

The distances between the test points are generated fully automatically to the minimum standard distance. If the distance is less than required by the standard, the affected test points are highlighted accordingly.

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Time-saving test mode 'Complete all indentations – then evaluate' and 'Optimum Stop' to complete test series as soon as the lower hardness limit has been undercut.

SIMPLIFIED LENS SELECTION


Based on the selected method (e.g. HV10), the suitable hardness range is displayed for each lens, which can be measured. The most suitable lens is also highlighted.

DUAL SYSTEM


With the Qpix Control 2 software, several QATM devices (for example a Qeye 800 and a Qness 60 A+ EVO) can be operated with the same PC system.

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

MODULE FOR SINTERED COMPONENTS

The software automatically searches for a position in the selected area where a hardness test indentation can be placed so that it is placed at a suitable location.

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



ANALISI STRUTTURALE RESA FACILE

MODULI SOFTWARE QPIX INSPECT

The intuitive and user-friendly Qpix INSPECT software functionality provides a comprehensive toolbox for microscopic evaluations and result documentation. The multifunctional software can be customized for customer-specific measuring tasks and complemented with add-on modules.

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INSPECT WELD SEAM MEASUREMENT

- | Standardised measurement and evaluation of weld seams
- | Automatic good/bad evaluation and report generation.

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IOT - INTERNET OF THINGS

THE PLATFORM FOR REMOTE ACCESS TO YOUR DEVICES

All QATM hardness testers with QpixControl2 and QpixT2 software seamlessly integrate into the Verder Scientific IoT platform, providing enhanced functionality and seamless connectivity.

- | **Real-time Monitoring:** Monitor your machinery in real time, from anywhere in the world. This data-driven approach empowers you to make informed decisions with ease.
- | **Live Notifications:** Be ahead of the curve with immediate alerts and updates. Real-time notifications ensure you stay informed about your equipment's performance.
- | **Effortless Backup:** Simplify your data protection. Whether you need to back up a single device or an entire fleet, our platform streamlines the process.
- | **Automatic & Free Software Updates:** Verder Scientific IoT ensures your customers' machines are consistently equipped with the latest software.

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© QATM – qatm.com / Il contenuto può essere soggetto a modifiche o correzioni